

FlowSpoon® FlowCup®

Materni's ready
to use ETO free
breastfeeding
aids

Innovative breastfeeding aids designed to support
the mother and child from the very first minute



Innovation in breastfeeding support

*Materni's breastfeeding aids are produced and packed in clean room Class 6 isolator conditions: CBM (Clean By Manufacturing).
Materni's proprietary simple but efficient feeding aids support breastfeeding from the very first moment.
The Materni single use FlowSpoon® & FlowCup® are suitable medical devices for Oral Rehydration Therapy (ORT).*



FlowSpoon®
Soft, flexible & smooth edges
MEPA 0012

Colostrum administration first

The soft and flexible FlowSpoon® allows safe collection and efficient administration of colostrum milk and Oral Rehydration Solution (ORS).
The spoon is individually packed and protected until point of use.

Safe and controlled cupping

Materni's FlowCup® is a unique, simple and affordable cupping device for the administration of expressed breast milk and ORS.
FlowCup® helps to avoid nipple-teat confusion and allows feeding at the natural pace of the baby, on top of being cheiloschisis friendly. No leakage or spillage. No milk-saliva mixing. Its simple design and intuitive use offer feeding flexibility to both professionals and parents. FlowCup® is individually packed & ready to use.

We offer 2 designs:

The **FlowCup® M1** was originally designed with a bridge to help control the milk flow in order to follow the baby's pace.

The **FlowCup® M2** was developed in close collaboration with healthcare professionals. Its design permits administration till the last drop while still offering a slow and controlled release.



FlowCup® M1
2 compartments
MEPE0111

FlowCup® M2
1 compartment
MEPE0113



FlowCup® and FlowSpoon® are registered trademarks. FlowCup® is a Materni patent.

Materni safe & easy FlowSpoon® and Materni FlowCup® standard product range

Article	Article code	Screw thread	Microbiology	Items/Box
Materni single use FlowCup® M1 (bridge) *	MEPE0111	EU	CBM	140
Materni single use FlowCup® M2 (no bridge) *	MEPE0113	EU	CBM	140
Materni single use FlowSpoon®	MEPA0012	n/a	CBM	300

* Materni bottles available separately.